


## Example: Decarburization

This module measures the depth or width of decarburization (for example when carbon is burnt from the surface of a steel during heat treating). As the decarburized surface layer can influence the properties of the finished product, it is important to be able to analyze it.


Compared with traditional manual techniques, this method can provide a significant improvement in laboratory productivity. The method is based on ASTM E-1077.

### PROCEDURE:


1. Load the image for the decarburization measurement.
2. Select the appropriate saved calibration scales from the drop down list at the bottom of the Window
3. Click on the Decarb icon  either from the Toolbar or from the Menu-bar
4. Choose the radio tab for
  - Complete Decarb
  - Partial Decarb
  - Total Decarb
5. Activate PARTIAL DECARB. Draw a line from the starting point to the end point of the partial decarb zone on the captured image.
6. Activate COMPLETE DECARB. Draw a line from the starting point to the end point of the complete decarb zone on the captured image. A blue line displays the width of the partial decarb and the red line represents the complete decarb width.
7. Click on TOTAL DECARB. The total of both the partial and complete decarb is displayed in the results window for field 1.
8. Repeat steps 5-7 to add more fields.
9. If you only require the complete decarb measurement, activate COMPLETE DECARB.
10. Draw one or more lines on the complete decarb layer.
11. To save the report, click on REPORT.
12. A sample information window will appear. Fill the required information. Click on "Save Data For Report Tool"
13. Click on "TO EXCEL" to see the result in Excel.

### EXAMPLE of Grain Size – DECARBURIZATION


File > Open (select image) or Capture > TWAIN Acquire (for live image)



2. Select Calibration scale at bottom of the Window



3. Click on the Decarb icon  either from the Toolbar or from the Menu-bar



4. Choose the radio tab for

- Complete Decarb
- Partial Decarb
- Total Decarb

Decarb Measurement ASTM E112/E1382-91

Trace Method: **Total Decarb**

☒ Partial Decarb ☐ Complete Decarb

Help For Drawing Lines

Select the Complete/Partial Option

Draw Lines to find out Width.

For Drawing Lines Click left button of mouse on starting point and drag it to the end point.

Result

Field :  Width:

**Total Decarb**

Field	Complete(Micron)	Partial(Micron)

5. Activate PARTIAL DECARB. Draw a line from the starting point to the end point of the partial decarb zone on the captured image.

Decarb Measurement ASTM E112/E1382-91

Trace Method: **Total Decarb**

☒ Partial Decarb ☐ Complete Decarb

Help For Drawing Lines

Select the Complete/Partial Option

Draw Lines to find out Width.

For Drawing Lines Click left button of mouse on starting point and drag it to the end point.


Result

Field :  Width:


**Total Decarb**

Field	Complete(Micron)	Partial(Micron)
Field 1		60.526

6. Activate COMPLETE DECARB. Draw a line from the starting point to the end point of the complete decarb zone on the captured image. A blue line displays the width of the partial decarb and the red line represents the complete decarb width.



7. Click on TOTAL DECARB. The total of both the partial and complete decarb is displayed in the results window for Field 1.



8. Repeat steps 5-7 to add more fields.
9. If you only require the complete decarb measurement, activate COMPLETE DECARB.
10. Draw one or more lines on the complete decarb layer.
11. To save the report, click on REPORT.
12. A sample information window will appear. Fill the required information. Click on "Save Data For Report Tool"
13. Click on "TO EXCEL" to see the result in Excel.

Grainsize1 [Compatibility]



File Home Insert Page Layout Formulas Data Review View

Cut Copy Paste Format Painter Clipboard Font Alignment Number

V14 fx

**MICROMEASUREMENT TEST REPORT**

NAME : APPLICATION :  
EVALUATION DATE : OPERATOR :  
SAMPLE INFO ID : MICROSCOPE OBJ :

Decarb ASTM E1077

Fields	Complete	Partial	Total
Field 1	552.632 Micron	78.947 Micron	631.579 Micron
Field 2	473.684 Micron	57.895 Micron	531.579 Micron
Field 3	739.474 Micron	31.579 Micron	771.053 Micron