

PACE[®] **TECHNOLOGIES**

www.metallographic.com

ZETA 2000 **Dispenser**

Automated / Manual Abrasive Dispenser

- 4 (24v) dedicated, programmable peristaltic pumps.
Setting adjustments for:
 - Frequency
 - Duration
 - Start time
 - Total time
- 7-inch touchscreen interface
- Compatible with NANO polishing machines as well as competitor polishing devices
- Controlled priming for each pump and dispensing lines for more immediate application of the abrasive / lubricant
- Flushing sequence for changing of abrasive/lubricant options
- Flow 1.2~90 ml/min
- Flex line for easy adjustment of dispensing nozzles
- 110/220 V compatible
- Direct line water inlet
- Easy to use

The ZETA 2000 is a 4-fluid automated / manual dispenser used for controlling and regulating the application of abrasives such as diamond, alumina, colloidal silica and lubricant solutions onto the polishing surface. The ZETA 2000 dispensing frequency, duration, starting time and total time can be programmed for manual or automated polishing. Priming the pumps lines allow for immediate application of the abrasive as required for consistent and reliable polishing. To prevent build-up or settling of abrasives in the pump and dispensing lines and pumps (especially for dispensing colloidal silica) an easy to use "Flush All Lines" command button has been added to the control panel. In addition, the ZETA 2000 has a hard connection to a water supply for more efficient and cleaning flushing of all lines.



ZETA-2000 ABRASIVE DISPENSER

PACE[®]
TECHNOLOGIES
www.metallographic.com

Dispenser Selection

<input checked="" type="checkbox"/>	1	<input checked="" type="checkbox"/>	2
<input type="checkbox"/>	3	<input type="checkbox"/>	4

Flush All
Lines

Prime
All Lines

Back to
Menu

Dispenser Settings

Frequency	Duration
0:30	0:02
Start Time	Total Time
0:00	4:00

START

STOP

Features

- Microprocessor controlled
- 7-inch touchscreen interface
- Digital LED readout
- 4 independent peristaltic pumps
- Interfaces with the FEMTO-S autopolishing heads for automated abrasive fluid dispensing
- Capable of dispensing NANOMETER alumina and SIAMAT colloidal silica and subsequently back flushing the lines to keep them clean and reuseable

Applications:

- Diamond**
- Alumina**
- Colloidal Silica**
- Lubricants**
- Water**

3601 E. 34th St Tucson, AZ 85713 USA
Phone +1-520-882-6598 FAX +1-520-882-6599
<https://www.metallographic.com>
email pace@metallographic.com



Machinery directive 2006/42/EC
Low Voltage Directive 2006/95/EC
Electromagnetic Compatibility directive 2004/108/EC

Technical Specifications:

Benefits:	Specification:
Pumps	4 (24V) dedicated, programmable peristaltic pumps. Setting adjustments for: -Frequency -Duration -Start time -Total time
Flow rates	1.2 - 90 ml/min
Controls	7-inch programmable touchscreen interface
Cleaning	Water flushing
Electrical	110V / 220V
Dimensions (WxDxH)	9.5 x 10.6 x 12-inch (240 x 270 x 300 mm)
Weight	12 lbs (5.5 kg)

Diamond Suspensions

Description	Monocrystalline Diamond (500 ml)	Polycrystalline Diamond (500 ml)
0.05 micron diamond	n.a.	GPC-0105-500
0.1 micron diamond	n.a.	GPC-0110-500
0.25 micron diamond	n.a.	GPC-0125-500
0.5 micron diamond	n.a.	GPC-0150-500
1 micron diamond	GMA-1001-500	GPC-1001-500
3 micron diamond	GMA-1003-500	GPC-1003-500
6 micron diamond	GMA-1006-500	GPC-1006-500
9 micron diamond	GMA-1009-500	GPC-1009-500
15 micron diamond	GMA-1015-500	GPC-1015-500
30 micron diamond	GMA-1030-500	GPC-1030-500
45 micron diamond	GMA-1045-500	GPC-1045-500

NANOMETER Polycrystalline Alumina Slurries

Description	0.05 micron Polycrystalline Alpha Alumina Slurry	0.05 micron Polycrystalline Acidic Alumina Slurr
Part Number	NA-1005-32	NA-1020-32
Slurry pH	4	10

Colloidal Silica Suspensions

Description	Part Number
SIAMAT colloidal silica (0.06 micron)	SI-100A-32
SIAMAT Blue colloidal silica (0.05 micron)	SI-100B-32
Non-crystallizing colloidal silica (0.02 micron)	SI2-100A-32
CMP alumina-colloidal silica (pH 10)	CMP-1005-32
CMP slurry (high viscosity)	CMP2-1005-32

PACE[®]
TECHNOLOGIES
www.metallographic.com



3601 E. 34th St Tucson, AZ 85713 USA
Phone +1-520-882-6598 FAX +1-520-882-6599
<https://www.metallographic.com>
email pace@metallographic.com