

FEMTO 1100S / FEMTO 1500S Autopolishing

8 and 10-inch Polishers

Benefits:

- Individual Force**
- Central Force**
- Variable speed head**
- Easy to use**



Automated Polishing

- FEMTO 1100S is a 1-6 Individual specimens polishing head
- FEMTO 1500S is an Individual and Central specimen polishing head
- 7-inch touchscreen interface
- 0-300 rpm Variable Speed Head
- 1 hp high torque dynamic servo motor
- Single or Double Wheel Polisher
- Retractable water rinse
- Manual or Auto operation
- Contra or Complementary Direction
- User friendly / Easy to Use

The FEMTO 1100S and FEMTO 1500 polishing heads can be installed at the factory or as a future upgrade for the NANO-1000S and NANO-2000S polishers. The FEMTO 1100S polishing head can polish 1 to 6 specimens individually and independent of the other specimens. This allows for easy removal and inspection of the specimen in-between each polishing step.


For flatter specimen preparation results the FEMTO 1500S offers a central force polishing fixture in addition to being able to polish specimens individually.

Individual Force:


Individual specimen force can be applied using the Quick Release Chuck system by screwing down the outside of the coupler so as to lock the specimen mounting plate into a rigid plane. This set-up allows for the individual holder to be removed and cleaned. This holder also eliminates the need to replane the samples if they are removed from the holder.

Central Force:

The individual pistons apply the polishing force to the landing pads and pushes down the holder through the spring loaded central male coupler. For central polishing a minimum of three samples locked into the holder is required; however, flatness is then fixed over all the mounted specimens.



Central Force



Individual Force



Individual vs. Central Force Automated Polishing

Individual Force: The force is applied to each sample individually. The main advantage for single or individual sample preparation is the elimination of the planar grinding steps. The primary disadvantage is that the specimen plane is not fixed in a rigid position during grinding and can result in loss of flatness across the specimen.


Central Force: For central force holders the specimens are locked into a rigid position or fixed plane so that the specimens cannot rock during sample preparation, especially the planar grinding or coarser grinding steps. The advantage is that the flatness across all the specimen mounts is better than individual force sample preparation. The main disadvantage is that 1-2 more planar grinding steps are required using coarser more damaging abrasives.

Individual Force Automated Polishing

Individual Force: Individual specimen force can be applied using the Quick Release Chuck system by screwing down the outside of the coupler so as to lock the specimen mounting plate into a rigid plane. This set-up allows for the individual holder to be removed and cleaned. This holder also eliminates the need to replane the samples if they are removed from the holder.




Individual Force Pistons in the off position




Individual Force Pistons in the on position - Force is applied to each specimen individually

Central Force Automated Polishing

Central Force: The individual pistons apply the polishing force to the landing pads and pushes down the holder through the spring loaded central male coupler. For central polishing a minimum of three samples locked into the holder is required; however, flatness is then fixed over all the mounted specimens.




Central Force in the off position



Central Force in the on position - Force is applied to the rigid specimen sample holder

ADVANTAGES

| FEMTO 1100S | FEMTO 1500S |
|---|---|
| Individual force | Central and Individual force |
| Reduced number of preparation steps and lower cost for consumables | Flatter specimens |
| Sample can be re-polished without having to planarize specimens in the holder | Specimens can be polished in both the contra and complementary head / base directions |



DISADVANTAGES

| FEMTO 1100S | FEMTO 1500S |
|---|---|
| Can result in loss of flatness across the specimen if heavy duty grinding is required | Minimum of 3-sample are required in the holder |
| It is recommended that the grinding and polishing steps only be done at matching head and base speeds with the rotation being in the same direction | Typically requires 1-3 additional coarser grinding steps to planarize all specimens |
| Fewer grinding steps | Additional grinding steps increase both consumables cost and total preparation time |
| Less grinding and subsurface damage | Coarser grinding can result in more surface and subsurface damage to the specimen |
| Do not need to planarize multiple samples | Need to planarize specimens again if they are removed from the fixed holder |




Technical Specifications:

| Specifications | FEMTO 1100S | FEMTO 1500S |
|---------------------------|---|---|
| Application of Force | Individual pistons (1-6 samples) | Individual pistons (1-6 samples) Central force (3-6 samples) |
| Specimen Holder | Fixed Individual Holder | Quick Release Chuck for Interchangeable Holders |
| Head Speed | 0-300 rpm | 0-300 rpm |
| Force | 1.1 - 13.5 lbs (5 - 60N) | 1.1 - 13.5 lbs (5 - 60N) |
| Power Supply | 110 / 220V | 110 / 220V |
| Motor | 1 hp dynamic torque servo motor | 1 hp dynamic torque servo motor |
| Weight | 60 lbs (28 kg) | 60 lbs (28 kg) |
| Dimensions (W x D x H) | 10 x 18 x 22-inch (254 x 457 x 560 mm) | 10 x 18 x 22-inch (254 x 457 x 560 mm) |
| Part No. | FEMTO-1100S | FEMTO-1500S |

| FEMTO 1100S and FEMTO 1500S Individual Specimen Holders | Part Number |
|--|--------------------|
| Individual Specimen holder for FEMTO 1100S (6 slots for 1 to 1.5-inch diameter rings) | SH-1100 |
| Individual Specimen holder for FEMTO 1500S with Quick Release Chuck (6 slots for 1 to 1.5-inch diameter rings) | SH-1150 |
| 1-inch diameter rings (6/set) 25 mm diameter rings (6/set) | SR-0100 SR-25mm |
| 1.25-inch diameter rings (6/set) 30 mm diameter rings (6/set) | SR-0125 SR-30mm |
| 1.5-inch diameter rings (6/set) 40 mm diameter rings (6/set) | SR-0150 SR-40mm |
| Individual Specimen holder for FEMTO 1100S (3 slots for 2-inch diameter specimens) | SH-1200 |
| Individual Specimen holder for FEMTO 1500S with Quick Release Chuck (3 slots for 2-inch diameter specimens) | SH-1250 |
| 2-inch diameter rings (3/set) 50 mm diameter rings (3/set) | SR-0200 SR-50mm |

| FEMTO 1500S Central Specimen Holders | Part Number |
|---|-----------------|
| 1-inch / 25 mm central specimen holder (3-6 samples) | QRC-SH100A-1500 |
| 1.25-inch / 30 mm central specimen holder (3-6 samples) | QRC-SH125A-1500 |
| 1.5-inch / 40 mm central specimen holder (3-6 samples) | QRC-SH150A-1500 |
| 2-inch / 50 mm central specimen holder (3 samples) | QRC-SH200A-1500 |
| Central Force Loading plate for FEMTO 1500S central holders | CF-LP-1500 |

| FEMTO 1500S | Part Number |
|---|-------------|
| Quick Release Chuck for dual central / individual specimens | QRC |
| Central holder male coupler for use with quick release chuck | QRC-CFC-V2 |
| Individual holder male coupler for use with quick release chuck | QRC-IFC-V2 |

 Machinery directive 2006/42/EC
 Low Voltage Directive 2006/95/EC
 Electromagnetic Compatibility directive 2004/108/EC

PACE[®]
TECHNOLOGIES
www.metallographic.com


PACE Technologies
 3601 E. 34th St
 Tucson, AZ 85713 USA
 Phone +1 520-882-6598
 FAX +1 520-882-6599
www.metallographic.com
 email pace@metallographic.com



Central specimen holder



Individual specimen holder



Quick Release Chuck



Central specimen male coupler adapter



Individual specimen male coupler adapter