




Example: Length measurement

Step 1:


File > Open (select image) or Capture > TWAIN Acquire (for live image)




Step 2. Select Calibration scale at bottom of the Window




Step 3-4. Select measurement by clicking on Measurement icon  or
Measurement toolbar and choosing measurement.



Step 5. Select Length Measurement, which opens up the following box



Step 6. Make measurement by Clicking and dragging from point A to point B on image (use box in Measurement screen for magnified image of cursor positioning for more accurate positioning)
Repeat for multiple measurements



Steps 7-9. Generating Reports.

Send TO EXCEL to create a report in Microsoft Excel.

Grainsize1 [Compatibility]

File Home Insert Page Layout Formulas Data Review View

Cut Copy Paste Format Painter Clipboard

Arial 10 A A

B I U

Font

Alignment

Wrap Text

Merge & Center

General

\$ % , .00 .00

Number



O8

MICROMEASUREMENT TEST REPORT

NAME : APPLICATION :

EVALUATION DATE : OPERATOR :

SAMPLE INFO ID : MICROSCOPE OBJ :

Sr.No	Description	Results (Micron)
L1	Length	363.625 Micron
L2	Length	151.127 Micron
L3	Length	273.697 Micron
L4	Length	207.294 Micron
L5	Length	317.244 Micron
L6	Length	477.543 Micron
L7	Length	225.227 Micron

Sheet1 Sheet2 Sheet3